



Dey-Engley Neutralizing HiCynth™ Broth

M1062

Dey-Engley Neutralizing HiCynth™ Broth is used in disinfectant testing where neutralization of the antiseptics and disinfectants is important for determining its bactericidal activity.

Composition**

Ingredients	Gms / Litre
HiCynth™ Peptone No.1	5.000
HiCynth™ Peptone No.5	2.500
Dextrose	10.000
Sodium thioglycollate	1.000
Sodium thiosulphate	6.000
Sodium bisulphite	2.500
Lecithin	7.000
Polysorbate 80	5.000
Bromocresol purple	0.020
Final pH (at 25°C)	7.6±0.2

**Formula adjusted, standardized to suit performance parameters

*Chemically defined peptones

Directions

Suspend 39.02 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense into tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle And Interpretation

Dey-Engley Neutralizing Broth is formulated as per the procedure described by Engley and Dey (1). Dey -Engley Neutralizing Broth is especially suited for environmental sampling where neutralization of the chemical is important to determine its bactericidal activity. Dey-Engley Neutralizing HiCynth™ Broth is prepared by completely replacing animal or vegetable peptones with chemically defined peptones to avoid BSE/TSE risks associated with animal peptone. A strongly bacteriostatic substance inhibits the growth and reproduction of bacteria without killing them. These bacteria hold the ability to cause infection under favourable conditions. Dey-Engley Neutralizing Broth Base and Dey-Engley Neutralizing Broth has the same formula but the former does not containing the neutralizing components. The Dey-Engley Neutralizing Broth neutralizes a broad spectrum of antiseptics and disinfectants including quaternary ammonium compounds, phenolics, iodine and chlorine preparations, mercurials, formaldehyde and glutaraldehyde. Dey-Engley Neutralizing Broth is used for the neutralization and testing of antiseptics and disinfectants according to the procedure of Engley and Dey (1).

HiCynth™ peptone No. 1 and HiCynth™ Peptone No.5 provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients. Dextrose is an energy source. The present formulation incorporate neutralizing substances for almost all the active products used as antiseptics and disinfectants. Sodium bisulfite neutralizes aldehydes; sodium thioglycollate neutralizes mercurials; sodium thiosulfate neutralizes iodine and chlorine (1); lecithin neutralizes quaternary ammonium compounds; and polysorbate 80, a non-ionic surface-active agent, neutralizes substituted phenolics (2-5). Bromocresol purple is an indicator for dextrose utilization. Due to the high concentration of lecithin in the broth medium, turbidity cannot be used to detect growth. Therefore, bromocresol purple and dextrose are added to the medium. Those organisms that ferment dextrose will turn the medium from purple to yellow. Growth of *Pseudomonas* species, which do not ferment dextrose, can be detected by the formation of a pellicle on the surface of the broth (1).

Neutralization Test

For testing disinfectants, prepare two sets of test tubes, one containing 9 ml Dey-Engley Neutralizing Broth (MCD1062) and other with 9 ml Dey-Engley Neutralizing Broth Base (MCD187). Add 1 ml of disinfectant under test. Mix well and allow it to stand for 15 minutes. Inoculate 0.1 ml of 1:100,000 dilution of overnight broth cultures and incubate at 37°C for 48 hours. Growth is indicated by a colour change from purple to yellow or pellicle formation. Growth in Neutralizing Broth and no growth in Neutralizing Broth Base indicate neutralization of disinfectant. To check bactericidal activity, both broth tubes are inoculated on D/E Neutralizing Agar (MCD186). Positive growth from negative tubes of Neutralizing Broth Base indicates bacteriostatic substance while negative growth indicates a bactericidal disinfectant. All positive tubes should show growth on Dey-Engley Neutralizing Agar. The control disinfectants used in test procedure are 2% chlorine, 2% formaldehyde, 1% glutaraldehyde, 2% iodine, 2% phenol, 1/750 quaternary ammonium compounds, 1/1000 mercurials etc.

Quality Control

Appearance

Light yellow to bluish grey homogeneous free flowing powder

Colour and Clarity of prepared medium

Purple coloured, opalescent solution in tubes with slight precipitate

Reaction

Reaction of 3.9% w/v aqueous solution at 25°C. pH : 7.6±0.2

pH

7.40-7.80

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 40-48 hours.

Organism	Inoculum (CFU)	Growth
Cultural Response		
<i>Bacillus subtilis</i> ATCC 6633	50-100	luxuriant
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant
<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	luxuriant
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant
<i>Staphylococcus aureus</i> ATCC 25923	50-100	luxuriant
<i>Escherichia coli</i> ATCC 8739	50-100	luxuriant
<i>Staphylococcus aureus</i> ATCC 6538	50-100	luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and use freshly prepared medium. Use before expiry date on the label.

Reference

- Engley and Dey, 1970. Chem. Spec. Manuf. Assoc. Proc., Mid-Year Meet., p. 100.

Revision :00 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.