



Baird Parker Agar (Agar medium O)

ME043

Baird Parker Agar is recommended for the isolation and enumeration of coagulase positive Staphylococci from food, pharmaceutical and other materials in accordance with European Pharmacopoeia.

Composition**

Ingredients	Gms / Litre
Pancreatic digest of casein	10.000
Beef extract	5.000
Yeast extract	1.000
Glycine	12.000
Sodium pyruvate	10.000
Lithium chloride	5.000
Agar	20.000
pH after sterilization (at 25°C)	6.8±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 63.0 grams in 950 ml purified /distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes or as per validated cycle. Cool to 50°C and add aseptically 50 ml concentrated Egg Yolk Emulsion (FD045) and 1 vial of sterile 1% Potassium Tellurite solution (FD052). Mix well before pouring into sterile Petri plates.

Warning: Lithium Chloride is harmful. Avoid all bodily contact and inhalation of vapours. On contact with skin wash with plenty of water immediately .

Principle And Interpretation

This medium is cited as Agar medium O in European Pharmacopoeia, 2008 (1) and is recommended for isolation and enumeration of coagulase positive *Staphylococcus aureus* . This medium was developed by Baird-Parker (2,3) from the Tellurite-glycine formulation of Zebovitz et al (4) for isolation of *S. aureus* from foods. *Staphylococcus* species are common contaminants in food, dairy, pharmaceutical and cosmetics related products (9). This medium is recommended for sterility checking of materials to detect *S. aureus* . Baird Parker medium was reported to be the best medium for selective detection of coagulase positive and entero-toxicogenic *Staphylococcus* (5). This medium was found to be less inhibitory to *S.aureus* than other media, at the same time being more selective (6,7). Subsequently it was officially adapted by the AOAC and European Pharmacopoeia (1,8).

Beef extract, yeast extract and pancreatic digest of casein provides essential mineral, vitamin and other growth requirements. Sodium pyruvate protects injured cells and helps recovery. Lithium chloride and potassium tellurite inhibit most of contaminating microflora except *S.aureus* . Glycine, pyruvate enhances growth of *Staphylococcus* . With the addition of egg yolk the medium becomes yellow and opaque. Glycine neutralizes aldehyde, while egg yolk neutralizes phenolic compounds, if any, in the test samples.

Proteolytic bacteria produce a clear zone around the colony in egg yolk containing media also known as Lecithinase reaction. A clear zone and grey-black colonies on this medium are diagnostic for coagulase positive Staphylococci. Upon further incubation, an opaque zone is developed around colonies, which can be due to lipolytic activity. Identity of *S. aureus* isolated on Baird-Parker Agar must be confirmed with a coagulase reaction and deoxyribonuclease test. The sterility of product is confirmed by absence of growth of *S. aureus* on this medium.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Please refer disclaimer Overleaf.

Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium

Basal medium: Yellow coloured clear to slightly opalescent gel After addition of Egg Yolk Emulsion and Tellurite Emulsion: Yellow coloured opaque gel forms in Petri plates.

Reaction

After sterilization, reaction of 6.3% w/v aqueous solution. pH : 6.8±0.2

pH

6.60-7.00

Cultural Response

Growth Promotion is carried out in accordance with EP. Cultural response was observed after an incubation at 35-37°C for 18-72 hours. Recovery rate is considered as 100% for bacteria growth on Soyabean Casein Digest Agar.

Cultural Response

ME043: Cultural characteristics observed after incubation at 35-37°C for 18-72 hours. Recovery rate is considered as 100% for bacteria growth on Soyabean Casein Digest Agar.

Organism	Inoculum (CFU)	Growth	Lot value (CFU)	Recovery	Colour of colony	Lecithinase
Growth Promoting						
<i>Staphylococcus aureus</i> ATCC 6538	50 -100	luxuriant	25 -100	>=50 %	grey-black shiny	Positive, opaque zone around the colony
Additional Microbiological testing						
<i>Staphylococcus aureus</i> ATCC 25923	50 -100	luxuriant	25 -100	>=50 %	grey-black shiny	Positive, opaque zone around the colony
<i>Proteus mirabilis</i> ATCC 25933	50 -100	good - luxuriant	50 -100	>=50	brown - black	Negative
<i>Micrococcus luteus</i> ATCC 10240	50 -100	poor - good	15 -40	30 -40 %	shades of brown-black (very small)	Negative
<i>Staphylococcus epidermidis</i> ATCC 12228	50 -100	poor - good	15 -40	30 -40 %	black	Negative
<i>Bacillus subtilis</i> ATCC 6633	50 -100	none - poor	0 -10	0 -10%	dark brown matt	Negative
<i>Escherichia coli</i> ATCC 8739	50 -100	none- poor	0 -10	0 -10%	large brown black	Negative
<i>Escherichia coli</i> ATCC 25922	50 -100	none- poor	0 -10	0 -10%	large brown black	Negative
<i>Escherichia coli</i> NCTC 9002	50 -100	none- poor	0 -10	0 -10%	large brown black	Negative

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

Reference

1. European Pharmacopoeia, 2008, European Department for the quality of Medicines.
2. Baird-Parker, A.C. 1962, J. Appl. Bact., 25: 12.
3. Baird-Parker, A.C. and Davenport, E., 1965, J. Appl. Bact., 28: 390.
4. Zebowitz, E., Evans J.B. & Niven C.F., (1955), J. Bact; 70:686.
5. Niskanean A and Aalto M, App. Env. Microbiol., 1978, 35:1233
6. Tardio and Baer, 1971, J. Assoc. Off. Anal. Chem., 54:728.
7. Baer, 1971, J. Assoc. Off. Anal. Chem., 54:732.
8. J. Assoc. off. Anal. Chem, 1971, 54:401.
9. FDA Bacteriological Analytical Manual, 2008, 18th ed., AOAC, Washington, DC

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.