



Medium 7. MacConkey Broth

MM083

Medium 7. MacConkey Broth is used for the selective identification of *E.coli* from pharmaceutical products in accordance with Indian Pharmacopoeia 2014.

Composition**

Ingredients	Gms / Litre
Gelatin Peptone	20.000
Lactose	10.000
Dehydrated ox-bile	5.000
Bromo cresol purple	0.010
pH after sterilization (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Pancreatic digest of gelatin

Directions

Suspend 35.01 grams in 1000 ml purified/ distilled water. Heat if necessary to dissolve the medium completely. Dispense into test tubes with inverted Durham tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes or as per validated cycle

Note: Directions specified are as per the concurrent edition of pharmacopoeia in force. Specified expiry period corresponds to this. User must ensure its compatibility with the latest edition.

Principle And Interpretation

MacConkey Broth is a modification of MacConkey Medium (1). Childs and Allen (2) demonstrated the inhibitory effect of neutral red and therefore substituted it by the less inhibitory bromocresol purple dye. BCP is more sensitive in recording pH variation in the medium. Indian pharmacopoeia (3) has recommended this medium for the presumptive identification of coliforms from water and other materials of sanitary importance.

Gelatin peptone provides essential growth nutrients. Lactose is the fermentable carbohydrate. Dehydrated ox-bile inhibits gram-positive organisms. Bromo cresol purple is the pH indicator in the medium. Lactose fermenting organisms gives yellow colour to the medium due to drop in pH as a result of lactose fermentation. Lactose non-fermenting organisms like *Salmonella* and *Shigella* do not alter the appearance of the medium. Gas produced gets entrapped in inverted Durhams tube.

Transfer homogenate in Casein Soyabean Digest Broth (MM011) containing 1 gm or 1 ml of the preparation to be examined to 100 ml MacConkey Broth Incubation is carried at 43°-45°C for 24-48 hours. For further isolation subculture on MacConkey Agar (MM081). Growth of red generally non-mucoid colonies, sometimes surrounded by a reddish precipitation zone, indicates presence of coliforms.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Purple coloured clear to slightly opalescent solution in tubes

pH

7.10-7.50

Cultural Response

Growth Promotion is carried out in accordance with the method of IP. For organisms not specified in pharmacopoeia, cultural response was observed after an incubation at 30-35°C for 18-48 hours.

Growth promoting properties

Clearly visible growth of microorganism comparable to that previously obtained with previously tested and approved lot of medium occurs at the specified temperature for not more than the shortest period of time specified inoculating ≤ 100 cfu (at 42-44°C for ≤ 24 hours).

Inhibitory properties

Please refer disclaimer Overleaf.

No growth of the test microorganism occurs for the specified temperature for not less than longest period of time specified inoculating ≥ 100 cfu(at 42-44°C for ≥ 48 hours).

Cultural Response

Organism	Inoculum (CFU)	Growth	Acid	Gas	Incubation temperature	Incubation period
Growth promoting						
<i>Escherichia coli</i> ATCC 8739	50 -100	luxuriant	positive reaction, yellow colour	positive reaction	42 -44 °C	≤ 24 hrs
Inhibitory						
<i>Staphylococcus aureus</i> ATCC 6538	$\geq 10^3$	inhibited			42 -44 °C	≥ 48 hrs
Cultural Response Additional Microbiological testing						
<i>Escherichia coli</i> ATCC 25922	50 -100	luxuriant	positive reaction, yellow colour	positive reaction	30 -35 °C	18 -24 hrs
<i>Escherichia coli</i> NCTC 9002	50 -100	luxuriant	positive reaction, yellow colour	positive reaction	30 -35 °C	18 -24 hrs
<i>Enterobacter aerogenes</i> ATCC 13048	50 -100	luxuriant	positive reaction, yellow colour	positive reaction	30 -35 °C	18 -24 hrs
<i>Salmonella Choleraesuis</i> ATCC 12011	50 -100	fair-good	negative reaction	negative reaction	30 -35 °C	18 -24 hrs
<i>Staphylococcus aureus</i> ATCC 25923	$\geq 10^3$	inhibited			30 -35 °C	≥ 48 hrs

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. MacConkey A. T., 1900, The Lancet, ii: 20.
2. Childs E. and Allen, 1953, J. Hyg: Camb. 51:468-477.
3. Indian Pharmacopoeia, 2014 Ministry of Health and Family Welfare, Govt. of India.

Revision : 02 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.