



## Medium 20.Pseudomonas Agar Medium for Detection of Fluorescein

MM120

Pseudomonas Agar for Detection of Fluorescein is used for the detection of fluorescein production by *Pseudomonas* species, and is in accordance with Indian Pharmacopoeia, 2007

### Composition\*\*

Ingredients	Gms / Litre
Pancreatic digest of casein	10.000
Peptic digest of animal tissue	10.000
Anhydrous dibasic potassium phosphate	1.500
Magnesium sulphate, 7H <sub>2</sub> O	1.500
Agar	15.000
pH after sterilization ( at 25°C)	7.2±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 37.23 grams of dehydrated medium in 1000 ml purified/distilled water, containing 10 ml glycerin. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

Pseudomonas Agar (For Fluorescein) is based on the formula described by King et al (1) and modified as in the Indian Pharmacopoeia (2) for the detection of fluorescein production a water soluble, chloroform insoluble fluorescent pigment by *Pseudomonas* species and in microbial limit tests(3). *Pseudomonas* are ubiquitous in environment and is common causative agent of burn, skin and nosocomial infections. They are also common contaminant of pharmaceutical and cosmetics related preparations. *Pseudomonas* strains are reported to produce phenazine pigments like Pyocyanin- blue green redox-active secondary metabolite pigment, pyorubin-rust brown pigment, -oxyphenazine- a breakdown product of Pyocyanin, pyoverdin-a water soluble yellow green pigments also known as fluorescein.

This medium enhances the elaboration of fluorescein by *Pseudomonas* and inhibits the pyocyanin formation. The fluorescein pigment diffuses from the colonies of *Pseudomonas* into the agar and shows yellow fluorescent colouration. Some *Pseudomonas* strains produce small amounts of pyocyanin resulting in a yellow-green colouration. Peptic digest of animal tissue and pancreatic digest of casein provides the essential nitrogenous nutrients, carbon, sulfur and trace elements for the growth of *Pseudomonas*. These nutrients are also conducive to the production of fluoroescien. Peptone and phosphorous in the medium enhances the production of pyoverdin/fluorescein pigment. Dipotassium phosphate buffers the medium while magnesium sulphate provides necessary cations for the activation of fluorescein production.

Salt concentration exceeding 2% affects pigment production. UV illumination may be bactericidal, so make sure that there is good growth before placing culture under UV light (3). The presence of Pseudomonas is confirmed by a positive oxidase test.

### Quality Control

#### Appearance

Cream to yellow homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.5% Agar gel

#### Colour and Clarity of prepared medium

Yellow coloured clear to slightly opalescent gel forms in Petri plates

#### pH

7.00-7.40

#### Growth Promotion Test

Growth Promotion is carried out in accordance with the harmonized method of IP. Cultural response was observed after an incubation at 33-37°C for not less than 3 days. Recovery rate is considered as 100% for bacteria growth on Soyabean Casein Digest Agar.

#### Cultural Response

MM120: Cultural characteristics observed after incubation at 33-37°C for 18-48 hours. Recovery rate is considered as 100% for bacteria growth on Soyabean Casein Digest Agar.

Organism	Inoculum (CFU)	Observed Lot value (CFU)	Recovery	Characteristic colonial morphology	Fluorescence in UV light	Growth
<b>Test for Pseudomonas aeruginosa</b>						
<i>Pseudomonas aeruginosa</i> ATCC 9027	50 -100	35 -100	≥70 %	Generally colourless to yellowish	positive	positive
<b>Additional Microbiological Testing</b>						
<i>Pseudomonas aeruginosa</i> ATCC 27853	50 -100	35 -100	≥70 %	Generally colourless to yellowish	positive	positive

#### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

#### Reference

1. King, Ward and Raney, 1954, J.Lab. Clin. Med., 44 : 301.
2. Indian Pharmacopoeia, 2007, Government of India, Ministry of Health and Family Welfare, Publications and Information Directorate (CSIR), New Delhi.
3. MacFaddin J., 1985, Media for Isolation-Cultivation-Identification and Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

Revision : 2 / 2015

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.