



Medium 4. Enterobacteria Enrichment Broth- Mossel

MM287

Enterobacteria Enrichment Broth, Mossel is used for selective enrichment of *Enterobacteriaceae* from pharmaceutical products in accordance with Indian Pharmacopoeia.

Composition**

Ingredients	Gms / Litre
Pancreatic digest of gelatin	10.000
Dextrose monohydrate	5.000
Dehydrated ox bile	20.000
Disodium hydrogen phosphate, dihydrate	8.000
Potassium dihydrogen phosphate	2.000
Brilliant green	0.015
pH after heating (at 25°C)	7.2±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 42.93 grams of dehydrated culture medium in 1000 ml purified/distilled water. Dispense 120 ml amounts in 250 ml flasks or 9 ml amounts in tubes. Stopper with cotton plugs or loose fitting caps. Heat in free flowing steam or boiling water (100°C) for 30 minutes. Avoid overheating of the medium. DO NOT AUTOCLAVE.

Principle And Interpretation

Enterobacteria Enrichment Broth, Mossel was formulated by Mossel et al and is recommended as an enrichment medium for bile tolerant Gram negative bacteria in bacteriological examination of non-sterile products and other products in food and animal feeds (1,2) and is recommended by Indian Pharmacopoeia(3). The formulation is in accordance with the harmonized methodology of USP/EP/BP/JP(5,6,7,11)

Pancreatic digest of gelatin and glucose monohydrate allows the growth of most of the members of *Enterobacteriaceae*. Brilliant green and ox-bile, purified are the inhibitory agents for gram-positive bacteria. Phosphates act as a good buffering agent and neutralizes acids produced by lactose fermenters that otherwise would adversely affect the growth of the organism. Lactose negative, anaerogenic lactose-positive or late lactose fermenting *Enterobacteriaceae* are often missed by the standard Coli-aerogenes test. To overcome this problem, lactose is replaced by glucose in the media. Phosphates form the buffering system of the medium. The cells damaged while drying or low pH are resuscitated in well-aerated Tryptone Soya Broth (MM011) for 2 hours at 25°C prior to enrichment in EE Broth. The resuscitation procedure is recommended for dried foods (8), animal feeds (9) and semi-preserved foods (10). EE Broth is an enrichment broth and a loopful of the enriched sample is then subcultured onto Violet Red Bile Glucose Agar (MM581) after an initial incubation at 30-35°C for 24 hours. Typical pink colonies from MM581 are subcultured for biochemical confirmation by oxidase and fermentation reactions (4). Decimal dilutions of the food homogenate are used if the expected counts are high or else initial suspension is used.

Quality Control

Appearance

Light yellow to greenish yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Emerald green coloured, clear solution without any precipitate

pH

7.00-7.40

Growth Promotion Test

Growth Promotion is carried out in accordance with the harmonized method of IP. Cultural response was observed after an incubation at 30-35°C for specified time.

Growth promoting properties

Clearly visible growth of microorganism comparable to that previously obtained with previously tested and approved lot of medium occurs at the specified temperature for not more than the shortest period of time specified inoculating ≤ 100 cfu (at 30-35°C for ≤ 24 hours).

Inhibitory properties

No growth of the test microorganism occurs for the specified temp for not less than longest period of time specified inoculating ≥ 100 cfu (at 30-35°C for ≥ 48 hours).

Cultural Response

MM287: Cultural characteristics observed after incubation at 30-35 °C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Acid	Incubation temperature	Incubation period
Growth Promoting					
<i>Escherichia coli</i> ATCC 8739	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	≤ 24 hrs
<i>Pseudomonas aeruginosa</i> ATCC 9027	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	≤ 24 hrs
Inhibitory					
<i>Staphylococcus aureus</i> ATCC 6538	$\geq 10^3$	luxuriant		30 -35 °C	≥ 48 hrs
Additional Microbiological testing					
<i>Escherichia coli</i> ATCC 25922	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	24 -48 hrs
<i>Escherichia coli</i> NCTC 9002	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	24 -48 hrs
<i>Pseudomonas aeruginosa</i> ATCC 27853	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	24 -48 hrs
<i>Enterobacter aerogenes</i> ATCC 14028	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	24 -48 hrs
<i>Proteus mirabilis</i> ATCC 25933	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	24 -48 hrs
<i>Salmonella Enteritidis</i> ATCC 13076	50 -100	luxuriant	positive reaction, yellow colour	30 -35 °C	24 -48 hrs
<i>Shigella boydii</i> ATCC 12030	50 -100	luxuriant	negative reaction	30 -35 °C	24 -48 hrs
<i>Staphylococcus aureus</i> ATCC 25923	$\geq 10^3$	inhibited		30 -35 °C	≥ 48 hrs

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

- Mossel D.A.A., Visser M. and Cornellisen A.M.R., 1963, J. Appl. Bact., 26(3):444.
- Van Schothurst M., et al, 1966, Vet Med., 13(3):273.
- Indian Pharmacopoeia, 2010, Govt. of India, Ministry of Health and Family Welfare, New Delhi, India.
- Mossel D.A.A. and Ratto M.A., 1970, Appl. Microbiol., 20:273.
- The United States Pharmacopoeia, 2011, The United States Pharmacopoeial Convention. Rockville, MD.
- British Pharmacopoeia, 2011, The Stationery office British Pharmacopoeia
- European Pharmacopoeia, 2011, European Dept. for the quality of Medicines.
- Mossel D.A.A. and Ratto M.A., 1970, Appl. Microbiol., 20:273.
- Mossel D.A.A. and Shennan J.L. and Clare V., 1973, J. Sci. Fd. Agric., 24 : 499.
- Mossel D.A.A., Ratto M.A., 1973, J. Fd. Technol., 8 : 97.

11.Japanese Pharmacopoeia, 2008.

Revision : 1 / 2011



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.