



Mucate HiVeg™ Broth

MV1226

Intended Use

Recommended for identification of enteropathogenic *Escherichia coli* and *Salmonella* species from milk and milk products on the basis of mucate utilization.

Composition**

Ingredients	Gms / Litre
Mucic Acid	10.000
HiVeg™ peptone	10.000
Bromo thymol blue	0.024
Final pH (at 25°C)	7.4±0.1

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 20.02 grams in 1000 ml distilled water. Dissolve mucic acid by slowly adding 5N sodium hydroxide and stirring. Dispense in 5 ml amounts in screw-capped tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 10 minutes.

Principle And Interpretation

Mucate HiVeg™ Broth is prepared by vegetable peptones which are free of BSE/TSE risks associated with animal peptones. It is the modification of Mucate Broth which is prepared based on the formula originally developed by Kauffman and Petersen (1) and recommended by APHA (2) for identification of enteropathogenic *Escherichia coli* from milk and milk products. This medium can also be used as an aid in differentiation of *Enterobacteriaceae* especially within *Salmonella* genus (3). Mucic acid is a saccharolactic acid or also called as tetrahydroxyadipic acid and act as a sole carbon source in the medium.

It is fermented by enteropathogenic *Escherichia coli*, *Salmonella* Paratyphi B and also, *Klebsiella pneumoniae* to produce acid which makes the medium yellow as the pH indicator is bromothymol blue (4). If the medium remains blue-green the organisms being tested does not utilize the mucate.

HiVeg™ peptone supplies the necessary nutrients to the organisms. Transfer a loopful of 24 hour Tryptone HiVeg™ Broth (M463) culture to Mucate HiVeg™ Broth. Include one negative control tube as well because un-inoculated tubes of Mucate Broth have been observed to turn blue on standing. Incubate at 48 ± 1 hour at 35°C. A negative test result is indicated by a blue or unchanged colour in this broth. 90 % of the strains are mucate positive.

Type of specimen

Dairy samples

Specimen Collection and Handling

For dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (2). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations :

Further biochemical identification is required for confirmation of species.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Blue coloured clear solution

Reaction

Reaction of 2.0% w/v aqueous solution at 25°C. pH : 7.4±0.1

pH

7.30-7.50

Cultural Response

MV1226: Cultural characteristics observed after an incubation at 35 - 37°C for 24 - 48 hours.

Organism	Inoculum(CFU)	Growth	Colour of medium
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	luxuriant	yellow
<i>Klebsiella pneumoniae</i> 13883 (00097*)	50-100	luxuriant	yellow
<i>Salmonella</i> Paratyphi B ATCC 8759	50-100	luxuriant	yellow

Key:- (*) Corresponding WDCM numbers

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

- 1.Kauffmann, F. and Petersen, A. 1956. Acta. Pathol. Microbiol. Scand., 38(6).
- 2.Wehr, H.M. and Frank, J.H. 2004. Standard Methods for the Examination of Dairy Products. 17 ed.
- 3.Ewing. 1986. Edwards and Ewings Identification of Enterobacteriaceae. 4 ed. N.Y: Elsevier Science Pub. Co., Inc.
- 4.MacFaddin, J. F. 1985. Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria vol. 1. Baltimore: Williams and Wilkins.
- 5.Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 03 / 2018

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.