



PHM010

Phyto Pspo Agar Base

Selective medium for the detection *Pseudomonas syringae* pv. *porri* on the seeds of leek

Composition **:

Ingredients	Grams/Litre
Peptone special	5.00
Sucrose	20.00
Di-potassium hydrogen phosphate	0.50
Magnesium sulphate anhydrous	0.13
Agar	20.00
Final pH (at 25°C)	7.5

**Formula adjusted standard to suit the performance parameter

Direction.

Suspend 45.63 grams in 1000 ml distilled water containing 1 gram of boric acid. Heat to boiling to dissolve the medium completely. Sterilize the medium by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50 °C and aseptically add the rehydrated contents of one vial of CNF supplement (PHS006) Mix well and pour into sterile Petri plates .

Principle and Interpretation

Bacterial blight of leek is caused by *Pseudomonas syringae* pv. *Porri*, a gram negative bacterium. The symptoms generally include Symptoms were water-soaked lesions at leaf tips, which eventually expanded down the length of the leaf and resulted in brown, elongated, stripe-like lesions with yellow margins. (1). *Pseudomonas syringae* pv. *Porri* is usually seed borne and therefore this medium is suitable for the isolation of this organism from the seeds of leek.

Peptone special provide nitrogenous compounds, carbon, vitamin B complex and trace ingredients. Sucrose serves as a source of carbohydrate. The inorganic phosphates in the medium serve as buffers.

Quality Control :

Appearance of Powder

Cream to yellow coloured homogeneous free flowing powder.

Gelling

Firm comparable with 2.0% Agar gel.

Colour and Clarity of prepared medium

Blue coloured opalescent gel forms in Petri plates

Reaction

Reaction of 4.56% w/v aqueous solution is pH 7.50 at 25°C.

PHM010**Phyto Pspo Agar Base****Cultural Response:**

Cultural characteristics observed after an incubation at 30-32°C for 5-6 days with added CNTB supplement-1 (PHS005).

Organism (ATCC)	Growth	Colour of the Colony
<i>Xanthomonas hortorum pv carotae</i>	Luxuriant	Light yellow cream , Light brown to Peach yellow, glistening, round and about 2-4 mm in diameter
Staphylococcus aureus	Inhibited	-

References:

1. Koike ST , Barak JD, Henderson DM, Gilbertson RL, Bacterial blight of leek: a new disease in California caused by *Pseudomonas syringae* ,Plant Disease 1999, 83:165-170,

Storage and Shelf-life :

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.