



PHS001

CNC Supplement

Selective supplement for the detection of *Pseudomonas syringae* on seeds of beans and also for the detection of *Ps. porri*, *Ps. pisi* and *Ps.tomato* on seeds of resp. leek, pea and tomato

Formula: per vial sufficient 1000ml medium

Cephalexin monohydrate	80 mg
Cycloheximide	100 mg

Rehydrate the contents of 1 vial with 2 ml of 0.2 N NaOH . Mix well and aseptically add to sterile molten Phyto Boric Acid Peptone Agar Base (PHM001). Mix well and pour into sterile Petri plates.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.