



Technical Data

PHS014

CNV Supplement

Selective supplement for the detection of *Pseudomonas savastanoi* pv. *phaseolicola*, *Pseudomonas syringae* and *Xanthomonas axonopodis* in seeds of beans.

Formula: (per vial, sufficient for 1000 ml medium)

| | |
|--------------------------|-------|
| Cephalexin monohydrate | 80 mg |
| Nystatin | 35 mg |
| Vancomycin Hydrochloride | 10 mg |

Directions:

Rehydrate the contents of 1 vial with 2 ml of 0.2N NaOH . Mix well and aseptically add to sterile molten Phyto Tyrosine Peptone Agar Base (MT Medium) (PHM018). Mix well and pour into sterile Petri plates.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.