
Sterile Distilled Water

Product Code: R041

Quality Control:

Appearance

Colourless clear liquid.

Reaction

Reaction of solution at 25°C.

Sterility Testing

Sterility of Distilled water was checked by incubating at 30-35°C for 14 days.

Results

Absence of any turbidity after 14 days at 30-35°C.

Storage and Shelf Life:

Store at 15 - 30°C.

Shelf life of the product is 18 months.

Use before expiry date given on the product label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.