



HiVeg™ Infusion No. 2

RM192V

Principle And Interpretation

HiVeg Infusion No. 2 is veg infusion recommended for use in place of Meat Infusion powder. It is highly nutritious ingredient which can be used for cultivation of wide variety of microorganisms. It can be used in standard nutrient media and as an additive in vaccine preparation in place of Meat Infusion powder.

Quality Control

Appearance

Light yellow to yellow, may have a slight greenish tinge Homogenous Free flowing powder, having Characteristic odour of protein, derived from vegetable source.

Solubility

Freely soluble in distilled/ purified water, insoluble in alcohol.

Clarity

1% w/v aqueous solution is clear without any haziness after autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Reaction

Reaction of 2% w/v aqueous solution at 25°C.

pH

5.50- 7.50

Microbial Load:

Total aerobic microbial count (cfu/gm)

By plate method when incubated at 30-35°C for not less than 3 days.

Bacterial Count : <= 2000 CFU/gram

Total Yeast and mould count (cfu/gm)

By plate method when incubated at 20-25°C for not less than 5 days.

Yeast & mould Count : <= 100 CFU/gram

Test for Pathogens

1. E.coli-Negative in 10 gms of sample
2. Salmonella species-Negative in 10 gms of sample
3. Pseudomonas aeruginosa-Negative in 10 gms of sample
4. Staphylococcus aureus- Negative in 10 gms of sample
5. C.albicans- Negative in 10 gms of sample
6. Clostridia- Negative in 10 gms of sample

Indole test

Tryptophan content: Passes

Cultural response

Cultural response observed after incubation at 35 - 37°C for 18-48 hours by preparing Nutrient HiVeg Agar (MV001), using HiVeg Infusion No.2 as an ingredient.

Cultural Response

Organism	Growth
Cultural response	
<i>Escherichia coli</i> ATCC 25922	Luxuriant
<i>Pseudomonas aeruginosa</i> ATCC 27853	Luxuriant
<i>Staphylococcus aureus</i> ATCC 25923	Luxuriant
<i>Salmonella Typhi</i> ATCC 6539	Luxuriant
<i>Streptococcus pyogenes</i> ATCC 19615	Luxuriant

Chemical Analysis

Total Nitrogen	$\geq 11.0\%$
Amino Nitrogen	$\geq 3.80\%$
Sodium chloride	$\leq 5.0\%$
Loss on drying	$\leq 7.0\%$
Residue on ignition	$\leq 20.0\%$

Storage and Shelf Life

Store below 30°C. Use before expiry date on the label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.