

RM612 Trypsin, From bovine pancreas**Product Number Packing**

RM612	: 1G
RM612	: 250MG


Product Information

Product Code	: RM612
Product Name	: Trypsin, From bovine pancreas
CAS No.	: 9002-07-7
EC No.	: 232-650-8
HS Code	: 3507 90 90
EC Index No.	: 647-010-00-7
Enzyme Classification No.	: 3.4.21.4
Shelf Life	: 3 years


Technical Specification

Appearance	: White to yellow crystals or powder
Solubility	: 33.3 mg soluble in 1 mL of water
Unit definition	: One BAEE unit will produce a ΔA_{253} of 0.001 per min at pH 7.6 at 25°C using BAEE as substrate. Reaction volume = 3.2 ml (1 cm light path)
Limit of Chymotrypsin	: ≤ 50 USP U/mg
Activity (on dry basis)	: ≥ 2500 USP U/mg

GHS Safety Information

Hazard Statement(s)	: H315-H319-H334-H335
Precautionary Statement(s)	: P261-P305+P351+P338-P342+P311
Signal Word	: Danger
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 36/37/38-42
S-Phrase(s)	: 22-24-26-36/37
WGK	: 1
RTECS	: YN5075000
Storage Temperature(°C)	: Store at 2 - 8°C
Hazard Symbol(s)	

Transport Information



Product Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein