



Product Information

RM6186 Protease (4000units/mg)

Product Number Packing

RM6186 : 100MG

Product Information

Product Code : RM6186
Product Name : Protease (4000units/mg)
Synonym : Protease E; Actinase E; Pronase E
CAS No. : 9036-06-0
EC No. : 232-909-5
Shelf Life : 3 years

Technical Specification

Appearance : White to tan powder.
Solubility : Dissolves in 0.01M Sodium acetate with 0.005 M Calcium acetate at pH 7.5.
Activity : ~ 4000 Units/ mg

Risk and Safety Information

S-Phrase(s) : 22-24/25
WGK : 2
RTECS : UK9595000
Storage Temperature(°C) : Store at 2 - 8°C

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein