

Horizontal Autoclave Al0052



REGISTERED OFFICE - 23, Vadhani Indl Est, LBS Marg, Mumbai - 400 086, India.
Tel: +91-22-6116 9797 / 2500 1607 | Fax: +91-22-2500 2286

CORPORATE OFFICE - A-516, Swastik Disha Business Park, Via Vadhani Indl Est, LBS Marg, Mumbai - 400 086, India.
Tel: +91-22-6147 1919 / 2500 3747 | Fax: +91-22-6147 1920 / 2500 5764 | Email: info@himedialabs.com
Web: www.himedialabs.com



... expect only quality from us™
CIN : U85195MH1982PTC028194



Specification	Value
Pressure Range	30 PSI
Temperature Range	135 Degree C
Power Consumption (KW)	2-6 KW
Chamber Volume (Litre)	200 Litre
Material	Stainless Steel
Capacity	1000 P.A Insulation
Wall	Double Wall Horizontal
Automation Grade	Fully Automatic
Size	Min 35 , Max 200 Ltr
Mount Type	Floor Mounted Autoclave
Approval Certificate	CE Certificate, ISO 9001 2008, WHO Certificate, ISO 9001 2008
Warranty	1 Year
Installation	Yes
Country of Origin	Made in India
Grade	SS304
Pressure inside the chamber	Variable from 5 psi. to 20 psi. pressure
Sterilization Process	1.2 kgf/cm sq.g (15 psi - 17 psi)
Working Temperature	121° C
VA106	550 x 750mm, 152 litre
Material	Stainless Steel
Automation Grade	Semi-Automatic, Fully Automatic
Temperature	121 degree
Power	220/230Volts, 50/60Hz AC supply
Operating Range	15psi to 20psi pressure

REGISTERED OFFICE - 23, Vadhani Indl Est, LBS Marg, Mumbai - 400 086, India.
Tel: +91-22-6116 9797 / 2500 1607 | Fax: +91-22-2500 2286

CORPORATE OFFICE - A-516, Swastik Disha Business Park, Via Vadhani Indl Est, LBS Marg, Mumbai - 400 086, India.
Tel: +91-22-6147 1919 / 2500 3747 | Fax: +91-22-6147 1920 / 2500 5764 | Email: info@himedialabs.com

Web: www.himedialabs.com



... expect only quality from us™
CIN : U85195MH1982PTC028194