

SEROLOGICAL WATER BATH



SEROLOGICAL WATER BATH

Serological water baths, with precise temperature control, are essential in labs and medical settings. They support microbiology, molecular biology, clinical labs, tissue culture, pharmaceutical research, and food safety testing. Their stability ensures reliable results across scientific disciplines.

SPECIFICATION

Inner side horizontal stainless-steel water bath with gabled lid

Easy Chamber Cleaning

Leak Free Inner Seamless Die Pressed Chamber.

Inner dimension W/D/H 500X300X150 mm

Safety alarm

Comes with Elegant looks

Full SS 304 Construction

Accurate PID Control with Dual Line Display and Soft Key

Illuminated Rocker Switch for On / Off

Temperature Range of Ambient +5 to 99.9 C°

High Stability of $\pm 0.3^{\circ}\text{C}$

Externally Cladder Heater

Better Heat Distribution

Minimum 2 years warranty

DOCUMENTS

