

# PET & PETG Bottles



**All-in-One Solution for Storage of Tissue Culture Media,  
Sera, Buffers, Reagents & Aqueous solutions**

**Choose Best. Choose Us.**

Because we value your samples just as much as you do

# 01

## STRINGENT INSPECTION AND TESTING

- ASTM standard D4991-07'
- Sterility Assurance Level (SAL) of 10<sup>-6</sup>

# 03

## INDIVIDUALLY PACKED

- For easier identification and easy opening linear tear.
- Applicable to product code: TCP269I, TCP270I, TCP271I and TCP272I.

# 05

## STACKABILITY

- Shrink-wrapped tray pack allows maximum storage and shipping space.

# 02

## LEAK-PROOF & LEACH-PROOF SYSTEM

- Free (ADCF) virgin resins

# 04

## VOLUME CAPACITY

- Four volume capacities available 125ml, 250ml, 500ml and 1000ml

