


HiPer® SaliSure Kit - ML254

- **Rapid Forensic Saliva Detection** – Results in 20-25 minutes using an enzymatic- chemical reaction.
- **Clear Visual Indication** – White spot appears on a purple background on detection paper.
- **High Sensitivity** – Detects saliva even at 1:1000 dilution and from month-old stains.
- **Versatile Surface Compatibility** – Detects Saliva from fabric, metal, plastic, glass, tiles, and wood

**Forensic Applications****On-site testing****Portable****Ready-to-use****Scan to view Video**

HiPer[®] SaliFind Kit - ML255



Positive Results
(With known Saliva)



Negative Results
(With distilled water)

Scan to
view Video



- Designed for rapid mapping of saliva traces for forensic purposes
- **Rapid Saliva Locator** – locate saliva traces within 20-25 minutes with a simple and rapid method
- **Clear Visual Indication** – White spot appears on a purple background on detection paper
- **High Sensitivity & Specificity** – Detect saliva diluted up to 1:1000 and identify stains up to one month old on various surfaces
- **Versatile Surface Compatibility** – Saliva can be mapped from fabric, metal, plastic, glass, tiles, and wood

Forensic Applications



On-site testing



Portable



Ready-to-use

