


PCR Plate

384 well Full-skirt White plate

PR29 - 384 well Full-skirt White

Feature : To reflect light & enhance the sensitivity of fluorescence signal analysis & quantification of amplified DNA.



PR30 - 384 well Full-skirt White plate with clear wells

Feature : For high sensitivity & precise optical measurements

Specifications (PR29 & PR30):

- ♦ Sterile PCR Plate
- ♦ Capacity : 40 μ L
- ♦ DNase/ RNase free
- ♦ Compatible HiMedia's Thermal cyclers : LA1111

