



Unlock Precision and Reliability with Hi-PCR® MatchFinder Kits

Multiplex Short Tandem Repeats (STR) Kits for Human DNA profiling for forensic analysis, relationship testing & research.



HiMedia's STR Portfolio

Feature	MBSTR001	MBSTR002	MBSTR003	MBSTR004
Name	Hi-PCR® Match Finder Autosomal 25ID Kit	Hi-PCR® Match Finder Autosomal 33ID Kit	Hi-PCR® Match Finder Y-43ID Kit	Hi-PCR® Match Finder Autosomal Mini25ID Kit
Target	Autosomal	Autosomal	Male	Autosomal
Number of Loci	25	33	43	25 in two panels
Dye System	6-dye	6-dye	6-dye	6-dye
Quality control markers	Yes	Yes	Yes	No
Sensitivity	Superior	Superior	Superior	Superior
Key Applications	Forensic, Genetic	Forensic, Genetic	Forensic, Genetic	Forensic, Genetic
Robustness to Inhibitors	Superior	Superior	Superior	Superior

Get the results you need, with the reliability you trust.