



Anoxomat® III Accessories

Anoxomat provides a stable environment for bacterial growth.

The system uses a disposable catalyst for anaerobic environments, and has a wide assortment of jars and petri dish holders.

Essential for Maximum Performance. Designed with you in mind.

Disposable catalyst saves time

Palladox™ is flexible and convenient so you work smarter when creating anaerobic environments.

Single-use, easy-to-use disposable sachet streamlines your workflow. You **don't waste precious time** reheating the catalyst and there is **no oven needed**.

- **Easy:** Made to fit and compatible with all jars
- **Saves time:** No reheating; cuts wait time by up to 2 hours
- **Simplifies workflow:** No guesswork on catalyst activity
- **Expedites processing:** So you optimize your time



Our accessories
provide an optimal
environment for
bacterial growth.

Choose from a wide assortment of jars

Connects to a variety of standard and ergonomic jars that hold different numbers of stacked culture plates of various sizes.



Maximize space with ergonomic jars

Designed for comfort and usability, you'll love the flexibility of these lightweight and compact jars.

- Comfort lid easily snaps on to create a secure seal
- Stack and nest to fit more jars in incubator
- Easy-grip handle makes transport a breeze
- Lid easily connects to gas supply
- Holds 12 plates or 24 plates per jar

Large capacity jars for large labs

Large laboratories need large capacity jars. Our standard jars do the trick. Large enough to hold up to 36 petri dishes, these jars meet the needs of large laboratories.

- Snap-shut gas supply coupling for easy connecting
- Lid clamp ensures a secure seal
- Variety of autoclavable Petri dish holders



Order today

Ergonomic Jars

AJ9049	Holds two stacks of 6 Petri dishes, 9-10 cm diameter
AJ9050	Holds two stacks of 12 Petri dishes, 9-10 cm diameter

Standard Large Capacity Jars

AJ9025	Holds three Microtiter plates, 13 x 9 cm Holds one stack of 12 Petri dishes, 15 cm diameter plates only
AJ9028	Holds three stacks of 12 Petri dishes, 9-10 cm diameter

Disposable Catalyst for Anaerobic Environments

AN3146	Palladox™ disposable sachet (box of 40) compatible with all jars
--------	--

Petri dish holders

PH 1040	Holds 12 Petri dishes, 9-10 cm diameter
PH 1050	Holds three stacks of 12 Petri dishes, 10 cm diameter
PH 1060	Holds 6 Petri dishes, 9-10 cm diameter
PH 1070	Holds 10 Petri dishes, 14.5 cm diameter
PH 1080	Holds three stacks of 12 Petri dishes, 6 cm diameter
PH 1090	Holds one stack of microtiter plates, 13 x 9 cm diameter
TH 0000	Tube holder



HIMEDIA®

For Life is Precious

HiMediaLaboratories™

www.himedialabs.com

Plot No. C-40, Road No. 21Y, Wagle Industrial Area, Thane (West) - 400604,

Tel : +91-22-6147 1919 / 6116 9797 / 6903 4800 | Fax : +91-22-6147 1920 | Email : info@himedialabs.com