

LA1025


Plate Inspection Port

Introducing



- Designed For Quality checking of prepared Plates.
- Inspection and Visibility made easy.
- Detect matter with in the medium.
- Low cost illuminator.
- Less Strain full and economical Light Weight & easy to carry.
- Remote operation for changing + - intensity of Light.
- Intensity 24900 LUX.
- Working Hours 10000.
- Plate inspection port is high performance LED Light, can be remote controlled by supporting controller, to control the light's on/off, and adjust the brightness.

CCT REMOTE CONTROL



Operating Procedure:

For first time use, hold on power Setup button until panel light blinks, then you can change desire colour or adjust brightness. If it is already use once, on need to hold power setup button again and can be use directly.

HIMEDIA®

HiMedia Laboratories Pvt. Ltd.

Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA. Tel.: 00-91-22-61471919 / 6116 9797 / 69034800
Fax : 00-91-22-61471920 Email : info@himedialabs.com Web : www.himedialabs.com