



Compressed Air Sampler (Diffuser) with LA881 HiAirflow Economical **LA1147**



HiMedia Compressed Air Sampler (diffuser) with LA881HiAirflow Economical to test the quality of compressed air and gases used in cleanrooms. A compressed air source is connected and the flow meter regulated to achieve the required flow rate.

Compressed Air Sampler (Diffuser) Attachment **LA1148**



Product Features

- Fully transportable.
- 90 mm standard Petri Plates.
- Sampling period timed to obtain the required total sample volume.

Certifications: IQ, OQ and PQ documentation available.

Delivery information: Self-contained unit, consisting of connections, airflow meter, pressure gauge, tap, and funnel head complete with plate housing.

HiMedia Laboratories Pvt. Ltd.

Plot No. C-40, Road No. 21Y, Wagle Industrial Estate, Thane (W) - 400604, Maharashtra, INDIA.
Tel.: +91 22 6116 9797/6147 1919 Fax: 00-91-22-6147 1920
Email : info@himedialabs.com,

HiMEDIA®
For Life is Precious