


Hi Touch Air Sampler



Introducing

Responsive Touchscreen

Wi-Fi Enabled


Software Excellence



LA1152



LA1151



Hi Touch Air Sampler

Technical specifications

Motor Type	Miniature Radial Blower, 6 VDC
Display Interface	4.3" TFT LCD with Resistive Touch Screen
Sampling Modes	Single Run, Multi Run
Sample Volume	Variable from 100 to 1000 liters with Increment of 100 liters and maximum volume of 2000 liters
Sample start Delay	Variable from 5 secs to 1800 seconds
Sample flow rate	100 liters/mins +/- 5%
Media Plate	Size : 90 mm
Sample Data Record	Not Editable Records up to 1,00,000
Input/Output	Wi-Fi
Audible Alarm	Electronic Buzzer
Material	LA1151 - Anodized Aluminium / LA1152 - ABS
Weight	LA1151 - 2.75 Kg/LA1152 - 2.1KG
Battery Pack	14.8 VDC, 6000 mAh Li-ion
AC/DC Charger Power Supply	Input : 100-300 VAC, 50 Hz, 14.8V-2A
Battery Charging and life	3 Hours to Charge Battery Battery Backup : Approximately 10 hours @ 100 liters/mins
Operating Range	0°C - 40°C 10-80% RH; Indoor Use
Operating system	Linux
Sampling Head Holes	219 (Hole size 1 mm)
Calibration Frequency	12 Months

Ordering Information

Product	Code
Hi Touch Air Sampler, MET 100	LA1151-1NO
Hi Touch Air Sampler, ABS 100	LA1152-1NO
Accessories for Hi Touch Air Sampler LA1151 & LA1152	
Compressed Air Assembly (CGA)	LA1151CGA - 1NO
Charger	LA1151A - 1NO
Battery	LA1151B - 1NO
Stainless Steel Drilled Head	LA1151D - 1NO
Anodized Aluminum Aspirating Head	LA1151E - 1NO
Disposable Drilled Head	LA1151F - 1NO
Carrying Case (Military Grade)	LA1151G - 1NO
Software with Software Installation Kit for Hi Touch Air Sampler	LA1151H - 1NO



HiMedia Laboratories Pvt. Ltd.

Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Estate, Thane, (West) 400604, Maharashtra, INDIA. Tel.: 00-91-22-61471919 / 6116 9797 / 69034800 Fax : 00-91-22-61471920 Email : info@himedialabs.com Web : www.himedialabs.com



For Website