

---

**Quantitative Filter Paper**  
**Ashless Filters ( Ashless Grades-Ash 0.007%)**  
**54 equivalent to grade 41**

**0540-1100-100C**

---

### Product Identifier

HS Code	: 4823 20 00
Storage	: Below 30C
Shelf life	: 4 Years

---

### Technical Specification

Description	: White coloured, circular, ashless quantitative filter paper. These standard ashless grades are high purity filter papers suitable for routine gravimetric techniques and in the preparation of samples for use in analysis involving instrumental techniques
Size	: 11 cm
Thickness	: 0.22 mm
Retention	: 10 um
Filtration Speed Rapidity	: 200 ml /min.
GSM	: 85 gm/m <sup>2</sup>
Loading capacity	: Medium
Use	: Used for coarse particulates and gelatinous precipitates such as iron, zirconium and aluminium hydroxides, and cobalt sulfide. also used in air pollution testing to determine gaseous compounds and fast analysis procedures for silicon

---

### Safety Information

UN No.	: Not dangerous good
Class	: -
Packing Group	: -