

---

**Quantitative Filter Paper**  
**Ashless Filters ( Ashless Grades-Ash 0.007%)**  
**54 equivalent to grade 41**

**0540-1250-100C**

---

**Product Identifier**

HS Code : 4823 20 00  
Storage : Below 30C  
Shelf life : 4 Years

---

**Technical Specification**

Discription : White coloured, circular, ashlessquantitative filter paper. These standardashless grades are high purity filter papers suitable for routine gravimetric technique sand in the preparation of samples for usein analysis involving instrumental techniques

Size : 12.5 cm  
Thickness : 0.22 mm  
Retention : 10 um  
Filtration Speed Rapidity : 200 ml/ min.  
GSM : 85 gm/m<sup>2</sup>  
Loading capacity : Medium

Use : Used for coarse particals and gelatinous precipitates such as iron, zirconium and aluminium hydroxides, and cobalt sulfide. also used in air pollution testing to determine gaseous compounds and fast analysis procedures for silicon

---

**Safety Information**

UN No. : Not dangerous good  
Class : -  
Packing Group : -