
Quantitative Filter Paper Ashless Filters (Ashless Grades-Ash 0.007%) 54 equivalent to grade 41

0540-1500-100C

Product Identifier

HS Code	: 4823 20 00
Storage	: Below 30C
Shelf life	: 4 Years

Technical Specification

Discription	: White coloured, circular, ashlessquantitative filter paper. These standardashless grades are high purity filter papers suitable for routine gravimetric techniques and in the preparation of samples for usein analysis involving instrumental techniques
Size	: 15 cm
Thickness	: 0.22 mm
Retention	: 10 um
Filtration Speed Rapidity	: 200 ml /min.
GSM	: 85 gm/m ²
Loading capacity	: Medium
Use	: Used for coarse particals and gelatinous precipitates such as iron, zirconium and aluminium hydroxides, and cobalt sulfide. also used in air pollution testing to determine gaseous compounds and fast analysis procedures for silicon

Safety Information

UN No.	: Not dangerous good
Class	: -
Packing Group	: -