

---

## Quantitative Filter Paper Ashless Filters ( Ashless Grades-Ash 0.007%) 74 equivalent to grade 40

**0740-4657-100S**

---

### Product Identifier

HS Code	: 4823 20 00
Storage	: Below 30C
Shelf life	: 4 Years

---

### Technical Specification

Description	: White coloured, circular, ashless quantitative filter paper. These standard ashless grades are high purity filter papers suitable for routine gravimetric techniques and in the preparation of samples for use in analysis involving instrumental techniques.
Size	: 46 x 57 cm
Thickness	: 0.21 mm
Retention	: 2.0 um
Filtration Speed Rapidity	: 30 ml/min.
GSM	: 95 gm/m <sup>2</sup>
Loading capacity	: Medium
Use	: Used in gravimetric analysis of components in heavy industries such as construction, mining of steel, filtration of trace elements such as silver chromate, lead sulfate, zinc and ammonium hydroxides and Soil Analysis.

---

### Safety Information

UN No.	: Not dangerous good
Class	: -
Packing Group	: -