
Quantitative Filter Paper Ashless Filters (Ashless Grades-Ash 0.007%) 94 equivalent to grade 42

0940-4657-100S

Product Identifier

HS Code	: 4823 20 00
Storage	: Below 30C
Shelf life	: 4 Years

Technical Specification

Description	: White coloured, circular, ashless quantitative filter paper. These standard ashless grades are high purity filter papers suitable for routine gravimetric techniques and in the preparation of samples for use in analysis involving instrumental techniques.
Size	: 46 x 57 cm
Thickness	: 0.2 mm
Retention	: 1.5 um
Filtration Speed Rapidity	: 4 ml/min.
GSM	: 100 gm/m ²
Loading capacity	: Medium
Use	: Used in critical gravimetric analysis of very fine particles such as barium and lead sulfates, stannic and nickel sulfide, calcium oxalate and calcium fluoride

Safety Information

UN No.	: Not dangerous good
Class	: -
Packing Group	: -