

---

**Qualitative Filter papers**  
**Standard grades(Non-wet-strengthened)**  
**601 equivalent to grade 1**

**6010-0900-100C**

---

**Product Identifier**

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

---

**Technical Specification**

Description	:	White coloured, circular, standard (non-wet strengthened) qualitative filter papers. These grades contain a high level of alpha-cellulose that provides high purity but weakens when wet.
Size	:	9 cm
Thickness	:	0.18 mm
Retention	:	11 um
Filtration Speed Rapidity	:	55 ml / min
GSM	:	87 gm/m <sup>2</sup>
Loading capacity	:	Medium
Use	:	Used in filtration of coarse precipitate, Soil analysis, Seed testing, in air pollution to collect dust laden air and in food industry to separate foodstuffs from associated liquid.

---

**Safety Information**

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-