
Qualitative Filter papers
Standard grades(Non-wet-strengthened)
601 equivalent to grade 1

6010-1250-100C

Product Identifier

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

Technical Specification

Description	:	White coloured, circular, standard (non-wet strengthened) qualitative filter papers. These grades contain a high level of alpha-cellulose that provides high purity but weakens when wet.
Size	:	12.5 cm
Thickness	:	0.18 mm
Retention	:	11 um
Filtration Speed Rapidity	:	55 ml / min
GSM	:	87 gm/m ²
Loading capacity	:	Medium
Use	:	Used in filtration of coarse precipitate, Soil analysis, Seed testing, in air pollution to collect dust laden air and in food industry to separate foodstuffs from associated liquid.

Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-