
Qualitative Filter Papers Standard grades (Non-wet-strengthened) 610 equivalent to grade 5

6100-1250-100C

Product Identifier

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

Technical Specification

Description	:	White coloured, circular, standard (non-wet strengthened) qualitative filter papers. These grades contain a high level of alpha-cellulose that provides high purity but weakens when wet.
Size	:	12.5 cm
Thickness	:	0.2 mm
Retention	:	1.5 μ m
Filtration Speed Rapidity	:	4 ml / minute
GSM	:	100 gm/m ²
Loading capacity	:	Medium
Use	:	Used for water analysis, biological product and fine precipitates such as cold-precipitated barium sulfate. In buchner funnels and to clarify cloudy suspensions.

Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-