
Qualitative Filter Papers Standard grades (Non-wet-strengthened) 631 equivalent to grade 4

6310-0900-100C

Product Identifier

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

Technical Specification

Description	:	White coloured, Circular, Standard (non-wet strengthened) qualitative filter papers. These grades contain a high level of alpha cellulose that provides high purity but weakens when wet.
Size	:	9 cm
Thickness	:	0.205 mm
Retention	:	10.0 um
Filtration Speed Rapidity	:	200 ml / minute
GSM	:	92 gm/m ²
Loading capacity	:	High
Use	:	Used for the filtration of coarse particles and gelatinous precipitates such as ferric aluminium and magnesium hydroxides. Used in air monitoring and gas filtration when high flow rate is needed.

Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-