

---

## Qualitative Filter Papers Standard grades (Non-wet-strengthened) 631 equivalent to grade 4

**6310-1250-100C**

---

### Product Identifier

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

---

### Technical Specification

Description	:	White coloured, Circular, Standard (non-wet strengthened) qualitative filter papers. These grades contain a high level of alpha cellulose that provides high purity but weakens when wet.
Size	:	12.5 cm
Thickness	:	0.205 mm
Retention	:	10.0 um
Filtration Speed Rapidity	:	200 ml / minute
GSM	:	92 gm/m <sup>2</sup>
Loading capacity	:	High
Use	:	Used for the filtration of coarse particles and gelatinous precipitates such as ferric aluminium and magnesium hydroxides. Used in air monitoring and gas filtration when high flow rate is needed.

---

### Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-