

---

## Qualitative Filter Papers Standard grades (Non-wet-strengthened) 642 equivalent to grade 2

**6420-1250-100C**

### Product Identifier

---

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

### Technical Specification

---

Description	:	White coloured, circular, standard (non-wet strengthened) qualitative filter papers. These grades contain a high level of alpha cellulose that provides high purity but weakens when wet.
Size	:	12.5 cm
Thickness	:	0.19 mm
Retention	:	2.0 um
Filtration Speed Rapidity	:	30 ml / min
GSM	:	97 gm/m <sup>2</sup>
Loading capacity	:	Medium
Use	:	Used in general laboratory application, agricultural and environmental, Laboratories for plant growth trials, soil testing and air analysis.

### Safety Information

---

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-