
Qualitative Filter Papers Wet-Strengthened Grades 961 equivalent to grade 113

9610-1500-100C

Product Identifier

HS Code : 4823 20 00
Storage : Below 30° C
Shelf life : 4 Years

Technical Specification

Description : White coloured, circular, standard (non-wet strengthened) qualitative filter papers. These grades contain a small amount of chemically stable, insoluble binding resin that strengthens the natural fibre to fibre bonds under wet condition

Size : 15 cm
Thickness : 0.42 mm
Retention : 40 um
Filtration Speed Rapidity : 360 ml / min
GSM : 125 gm/m²
Loading capacity : Very high

Use : Used for the filtration of coarse and gelatinous precipitates

Safety Information

UN No. : Not dangerous good
Class : -
Packing Group : -