
Qualitative Filter Papers Wet-Strengthened Grades 961 equivalent to grade 113

9610-4657-100S

Product Identifier

HS Code	:	4823 20 00
Storage	:	Below 30° C
Shelf life	:	4 Years

Technical Specification

Description	:	White coloured, circular, standard (non-wet strengthened) qualitative filter papers. These grades contain a small amount of chemically stable, insoluble binding resin that strengthens the natural fibre to fibre bonds under wet condition
Size	:	46 x 57 cm
Thickness	:	0.42 mm
Retention	:	40 µm
Filtration Speed Rapidity	:	360 ml / min
GSM	:	125 gm/m ²
Loading capacity	:	Very high
Use	:	Used for the filtration of coarse and gelatinous precipitates

Safety Information

UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-