



Minimum Essential Medium Eagle (MEM) 2X

**With Earle's salts, NEAA, L-Glutamine and Sodium bicarbonate
Without Phenol red**

Product Code: AL178A

Product Description :

Minimum Essential Medium (MEM) is a modification of Basal Medium Eagle (BME). It was developed by Harry Eagle to meet the specific nutritional requirements of certain subtypes of HeLa cells and normal mammalian fibroblasts. MEM includes higher concentration of amino acids so as to closely approximate the protein composition of cultured mammalian cells. MEM can be used either with Earle's salts or Hanks' salts and can also be additionally supplemented with Non-essential amino acids (NEAA). This medium can be further modified by eliminating calcium to facilitate growth of cells in suspension cultures.

AL178A is 2X Minimum Essential Medium Eagle with Earle's salts, NEAA, L-glutamine and sodium bicarbonate. It does not contain phenol red. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition:

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	530.000
Magnesium sulphate anhydrous	195.440
Potassium chloride	800.000
Sodium bicarbonate	4400.000
Sodium chloride	13600.000
Sodium dihydrogen phosphate anhydrous	244.000
AMINO ACIDS	
Glycine	15.000
L-Alanine	17.800
L-Arginine hydrochloride	252.000
L-Asparagine monohydrate	30.000
L-Aspartic acid	26.600
L-Cystine dihydrochloride	62.600
L-Glutamic acid	29.400
L-Glutamine	584.000
L-Histidine hydrochloride monohydrate	84.000

L-Isoleucine	104.000
L-Leucine	104.000
L-Lysine hydrochloride	145.000
L-Methionine	30.000
L-Phenylalanine	64.000
L-Proline	23.000
L-Serine	21.000
L-Threonine	96.000
L-Tryptophan	20.000
L-Tyrosine disodium salt	103.800
L-Valine	92.000
VITAMINS	
Choline chloride	2.000
D-Ca-Pantothenate	2.000
Folic acid	2.000
Nicotinamide	2.000
Pyridoxal hydrochloride	2.000
Riboflavin	0.200
Thiamine hydrochloride	2.000
i-Inositol	4.000
OTHERS	
D-Glucose	2000.000

Quality Control:

Appearance

Clear solution, with yellow tinge.

pH

7.00 - 7.60

Osmolality in mOsm/Kg H₂O

575 - 615

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts and comparing it with a control medium.

Endotoxin Content

NMT 1EU/ml

Storage and Shelf Life:

Store at 2-8°C away from bright light.

Shelf life is 12 months.

Use before expiry date given on the product label.

Disclaimer :

Revision : 06/2025

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.