

BG001 Glass Slides, Ground Edges, Plain**Product Number****Packing**

BG001	: 10 x 50 No.
BG001	: 5 x 50 No.

Product Information

Product Code	: BG001
Product Name	: Glass Slides, Ground Edges, Plain

Technical Specification

Description	: Glass slides having plain ground edges, with 90°corners. The slides are of premium quality.
Dimensions	: 76 mm x 26 mm
Thickness	: ~0.95 mm

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods