



Product Information

CG147

Slide Staining Jar, With ground-in Glass Lid

Product Number**Packing**

CG147 : 5 x 1 No

Product Information

Product Code : CG147

Product Name : Slide Staining Jar, With ground-in Glass Lid

Technical Specification

Description : Slides Staining Jar with ground in glass lid for 30 pieces 25.4 x 76.2 mm (1" x 3") slides. It is available in transparent and amber coloured and suitable for stainless steel rack of CG147.

Use : Used in staining of slides.

Transport Information

ADR/RID : Not Dangerous Goods

IMDG : Not Dangerous Goods

IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein