

CG282E : PCR Tubes, Flat lid, sterile, Vol. 0.2ml

Product Number : Packing

CG282E : 1X1000NO

Product Information

Product Code : CG282E

Product Name : PCR Tubes, Flat lid, sterile, Vol. 0.2ml

Technical Specification

Description : Sterile/ Autoclavable PCR tubes having conical base with flat lid frosted for easy labelling

Clarity : Clear without any particles and fibers.

Color : Clear, natural

Material : Polypropylene

Graduations : No

Height of Tube : 19.5 - 20.5 mm

Height of Cap : 3.5 - 4.5 mm

Inner mouth Diameter of tube : 5.0 - 6.0 mm

Outer mouth Diameter of tube : 6.5 - 7.5 mm

Outer mouth Diameter of cap : 5.0 - 6.0 mm

Capacity : 0.2 mL

Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

Product alignment : Lid closure of the tube is tight.

Use : Suitable for PCR applications.


Please refer disclaimer Overleaf.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

Technical Assistance

At HiMedia we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at mb@himedialabs.com.



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Private Limited,
Reg. Off: Plot No. C-40, Road No. 21Y,
MIDC, Wagle Industrial Area, Thane,
(West) 400604, Maharashtra, INDIA.
Web: www.himedialabs.com



05/2025

PICG282E_1/0522

CG282E-01

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. Reg. office : Plot No. C-40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane, (West) 400604, Maharashtra, INDIA.
Customer Care No.: 00-91-22-6116 9797 Tel: 00-91-22-6147 1919, 6903 4800 Email: techhelp@himedialabs.com Website: www.himedialabs.com