
Pyrazinamide

CMS10335

Product Identifier

CAS No.	:	98-96-4
EC No.	:	202-717-6
Molecular Formula	:	C ₅ H ₅ N ₃ O
Molecular Weight	:	123.11
HS Code	:	2933 99 00
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to Off-White crystals or powder needles or chunks
Solubility	:	10 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	187 - 192°C
Assay (GC/Total nitrogen)	:	>= 98.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UQ2275000
WGK	:	3