

Actinomycin D, From Streptomyces parvalus, For Biochemistry

CMS1243


Product Identifier

CAS No.	:	50-76-0
EC No.	:	200-063-6
Molecular Formula	:	C ₆₂ H ₈₆ N ₁₂ O ₁₆
Molecular Weight	:	1255.42
Synonym	:	Dactinomycin
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Orange to red crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of ethanol
Assay (HPLC)	:	min. 90.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H300- H350- H360FD
Precautionary Statement(s)	:	P201- P280- P301 + 310- P308 + P313- P330- P501
UN No.	:	3462
Class	:	6.1
Packing Group	:	II
RTECS	:	AU1575000
WGK	:	3