

CMS1334 Ceftriaxone sodium salt**Product Number Packing**

CMS1334 : 1G

Product Information

Product Code : CMS1334
Product Name : Ceftriaxone sodium salt
Synonym : Ceftriaxone disodium salt hemi(heptahydrate)
Molecular Formula : $C_{18}H_{16}N_8Na_2O_7S_3 \cdot 3.5H_2O$
Molecular Weight : 661.60
CAS No. : 104376-79-6
EC No. : 277-930-0
HS Code : 2941 90 00
Shelf Life : 3 years

Technical Specification

Appearance : White to light yellow crystals or powder
Solubility : 100 mg soluble in 1 mL of water
Water (K.F.) : 8.0 - 11.0%

GHS Safety Information

Hazard Statement(s) : H315-H317-H319-H334-H335
Precautionary Statement(s) : P261-P280-P305+P351+P338-P342+P311
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38-42/43
S-Phrase(s) : 26-36
WGK : 1
Storage Temperature(°C) : Store at 2 - 8°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods