

---

## p-Aminobenzoic acid, Hi-LR<sup>®</sup>

**CMS171**

### Product Identifier

---

CAS No.	:	150-13-0
EC No.	:	205-753-0
Molecular Formula	:	C <sub>7</sub> H <sub>7</sub> NO <sub>2</sub>
Molecular Weight	:	137.14
Synonym	:	PABA, Hi-LR <sup>®</sup> ; Vitamin B <sub>10</sub> , Hi-LR <sup>®</sup>
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	3 years

### Technical Specification

---

Appearance	:	White to beige crystals or powder or chunks or solid
Solubility	:	33.3 mg soluble in 1 mL of ethanol
FTIR	:	Matches with the standard pattern
Melting range	:	185 - 188°C
Assay (NaOH Titration/GC)	:	99.00 - 102.00%

### Safety Information

---

Hazard Statement(s)	:	H412
Precautionary Statement(s)	:	P273- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	2