

Retinol acetate

CMS183


Product Identifier

| | | |
|-------------------|---|--|
| CAS No. | : | 127-47-9 |
| EC No. | : | 204-844-2 |
| Molecular Formula | : | C ₂₂ H ₃₂ O ₂ |
| Molecular Weight | : | 328.49 |
| Synonym | : | Vitamin A acetate powder |
| HS Code | : | 2936 21 00 |
| Storage | : | 2 - 8°C |
| Shelf life | : | 3 years |

Technical Specification

| | | |
|------------|---|------------------------------------|
| Appearance | : | White to yellow crystals or powder |
| Solubility | : | Insoluble in water |
| Potency | : | 475,000 to 575,000 IU/gm |

Safety Information

| | | |
|----------------------------|---|--|
| Hazard Pictogram(s) | : |  |
| Signal Word | : | Danger |
| Hazard Statement(s) | : | H315- H361 |
| Precautionary Statement(s) | : | P201- P280- P302 + P352- P308 + P313- P332 + P313- P501 |
| UN No. | : | Not available |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | VH6825000 |
| WGK | : | 2 |