

## Vitamin E

CMS185


### Product Identifier

CAS No.	:	10191-41-0
EC No.	:	233-466-0
Molecular Formula	:	C <sub>29</sub> H <sub>50</sub> O <sub>2</sub>
Molecular Weight	:	430.71
Synonym	:	(±)-α-Tocopherol, liquid
HS Code	:	2936 28 00
Storage	:	2 - 8°C
Shelf life	:	3 years

### Technical Specification

Appearance	:	Colourless to brown viscous oily liquid
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Assay (HPLC)	:	>= 98.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H317
Precautionary Statement(s)	:	P261- P280- P321- P333 + P313- P362 + P364- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	GA8746000
WGK	:	1