
Vitamin E acetate

CMS1870

Product Identifier

CAS No.	:	7695-91-2
EC No.	:	231-710-0
Molecular Formula	:	C ₃₁ H ₅₂ O ₃
Molecular Weight	:	472.74
Synonym	:	DL- α -Tocopherol acetate
HS Code	:	2936 28 00
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Colourless to greenish-yellow viscous liquid
Solubility	:	100 mg soluble in 1 mL of chloroform
Refractive index (n _{20/D})	:	1.4940 - 1.4990
Density (at 25°C)	:	0.910 - 1.010 g/mL
Assay (GC/HPLC)	:	min. 96.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1