

Ethambutol dihydrochloride

CMS1890


Product Identifier

CAS No.	:	1070-11-7
EC No.	:	213-970-7
Molecular Formula	:	C ₁₀ H ₂₄ N ₂ O ₂ ·2HCl
Molecular Weight	:	277.23
HS Code	:	2922 19 90
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to pale yellow crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
Melting range	:	~202°C
Assay (NT/AT/HPLC)	:	99.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H360
Precautionary Statement(s)	:	P201- P280- P308 + P313- P405- P501
UN No.	:	Not dangerous good
Class	:	-
Packing Group	:	-
RTECS	:	EL3854000
WGK	:	1