

Lincomycin hydrochloride

CMS211


Product Identifier

CAS No.	:	859-18-7
EC No.	:	212-726-7
Molecular Formula	:	C ₁₈ H ₃₄ N ₂ O ₆ S.HCl
Molecular Weight	:	443.0
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to pale yellow crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
Water (K.F.)	:	3.0 - 6.0%
Assay (HPLC/AT, anhydrous basis)	:	min. 95.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RH6315000
WGK	:	1

Other Information

1 million units per vial