

---

## Nystatin

**CMS212**

---

### Product Identifier

CAS No.	:	1400-61-9
EC No.	:	215-749-0
Molecular Formula	:	C <sub>47</sub> H <sub>75</sub> NO <sub>17</sub>
Molecular Weight	:	926.09
Synonym	:	Mycostatin
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

---

### Technical Specification

Appearance	:	Off-white to slightly brownish hygroscopic powder
Solubility	:	100 mg soluble in 1 mL of dimethylformamide
Potency	:	>= 4400 USP units/mg

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RF5950000
WGK	:	1

---

### Other Information

1,000,000 units per vial